

FastPak® TableTop

Oral Solid Packaging and Barcoding System



Compact unit dose barcoding/packaging solution for any size pharmacy requiring improved efficiency and safety.

An RxWorks™ Suite Solution

The FastPak® TableTop is one of many RxWorks Suite technology and consulting solutions that optimize the medication management process and enable improved patient care, increased cost savings, and market growth opportunities.

Flexible to meet the needs of any health system pharmacy

The FastPak® TableTop can be your primary unit dose packaging automation system, a complementary addition to other unit dose packaging systems or a readily available backup solution. It is well-suited to the medication packaging needs of any type of facility, especially hospitals and long-term care facilities. The bottom line is that your pharmacy is incomplete without a FastPak TableTop.

Feature-rich to ensure packaging safety and accuracy

The FastPak Tabletop software incorporates barcode scanning of medication bulk bottles and images of medications being packaged to ensure the accuracy of unit dose packaging.

Easy to learn, easier to operate

Your pharmacy will enjoy the efficiency of the FastPak TableTop within minutes of installation. The solution combines a touch screen monitor, on-board help videos and easy-to-learn, proprietary software for simple operation.

Reduces medication costs

Purchasing oral solid medications in bulk saves the typical hospital 4-6 cents per oral solid dose. These savings can help achieve a payback within one year in a typical 100-bed hospital, when used as a primary packaging system.

Packaging allows extended dating, color options

The FastPak TableTop boasts the smallest packaging on the market while still accommodating large tablets and/or gel caps. The size of our packaging also increases the capacity of your unit-based cabinets.


AmerisourceBergen®

The Best Medicine for Healthcare™

FastPak® TableTop

Oral Solid Packaging and Barcoding System



Enables bedside medication verification at the point of care

Barcode Medication Administration (BCMA) systems are only effective if the medications are barcoded when they arrive at the bedside. When you combine the FastPak TableTop and the Barcoding Station (optional) for barcoding injectable and unit of use items, the result is a complete barcoding solution. To further enhance an effective BCMA system, 2D barcodes can display expiration and lot number.

Additional Product Highlights:

- The FastPak TableTop oral solid medication packaging system is a state-of-the-art, affordable system for prepackaging medications into unit dose packages.
- Horizontal operation virtually eliminates pill crushing and gel cap melting that is common with vertical feed systems. The AutoPac optional autofeeding device is recommended for runs of 600+ packages.
- Its small footprint allows the system to be placed in any location with no disruption.
- The barcode feature enables pharmacies to implement bedside barcode verification for medication administration. Available barcode formats are 2D and linear.
- In addition to barcoding verification after the packaging process, you have the ability to barcode verify during the process to ensure readability before completion of large runs.

Specifications:

FastPak TableTop

Speed: Up to 80 barcoded packages per minute

Weight: 175 lbs. (optional AutoPac 30 lbs.)

Footprint: 4.5 sq. ft.

Dimensions: 14" D x 46" W x 24" H

Components:

FastPak TableTop software
Windows® computer with keyboard, mouse
Touch screen monitor
Barcode scanner
UPS for computer

Packaging: Poly in five different light-resistant colors, including clear

Options: AutoPac autofeeding device

Power Requirements: 110V, 50/60Hz

Barcoding Station (Optional)

AmerisourceBergen's® Barcoding Station is a cost-effective system that accurately barcodes all items not barcoded by the manufacturer. When the FastPak TableTop is combined with the Barcoding Station, any size hospital can achieve a complete barcode solution.

Components:

AMS
Windows® client computer with modem, keyboard, mouse
Touch screen monitor
Cordless barcode scanner
Barcode label printer
UPS for client computer

Power Requirements: 110V, 60Hz, 15 amps

For more about these and other RxWorks Suite solutions, contact AmerisourceBergen at 877-892-1254; email at solutions@amerisourcebergen.com or visit our website at www.AmerisourceBergenDrug.com.

AmerisourceBergen
1300 Morris Drive
Chesterbrook, PA 19087
877-892-1254

AmerisourceBergenDrug.com


AmerisourceBergen®

The Best Medicine for Healthcare™